

SODIUM HYALURONATE OVDs

ALSAVISC

Sodium Hyaluronate Viscoelastics



ALSAVISC

OVDs made for all
needs in cataract
surgery



Dispersive and Cohesive Biofermented sodium hyaluronate OVDs

Cohesive range / 1% - 1.4% - 1.6% - 1.8%

- Space creating
- Stable anterior chamber
- Easy implantation

Dispersive range / 3%

- Endothelial protection
- Controlled capsulorhexis
- Easy surgery




1

- Sodium hyaluronate from bacterial fermentation
- pH 6.8 – 7.4
- Osmolality 270-400 mOsm/kg

2

- 1 ml syringe
- Storage 2°C - 25°C



SODIUM HYALURONATE OVDs

ALSWOBR-001

ALSANZA

Cataract environment



VISCOELASTICS



KNIVES



IRRIGATION



PHACO



IOLs

Address:

Hermann-Burkhardt-Straße 3
72793 Pfullingen

Germany

Tel: +49 (0) 7121-69 065-20

Email: info@alsanza.com

www.alsanza.com

ALSANZA
Medizintechnik und Pharma GmbH